#### GC Technical Note

GT144 GL Sciences Inc.

## Volatile Component Analysis of Combustion Using The Portable Thermal Resolver HandyTD TD265

Assuming a fire site where fuels such as gasoline, kerosene, and light oil were used, we used HandyTD TD265 to perform a concentration analysis of the components that volatilize from the recovered material after combustion. In the recovered material, components derived from the fuel used for combustion were detected.

#### **Standard Solution Analysis**



#### <Standard Sample Preparation Procedure>

- 1. Dilute gasoline, kerosene, and light oil with hexane.
- 2. Add the diluted sample to the adsorption tube for thermal desorption and analyze by TD/GC/MS.

#### Condition

| <u>GC Conditions</u> |                                              |
|----------------------|----------------------------------------------|
| System               | : GC/MS                                      |
| Column               | : InertCap 1                                 |
|                      | 0.25 mm I.D. $\times$ 60 m df = 0.25 $\mu$ m |
| Col.Temp.            | : 40 °C (5 min) - 10 °C /min - 300 °C        |
| Carrier Gas          | : He 150 kPa                                 |
| Injection            | : Split 1:100                                |
|                      | 240 °C                                       |
| Detection            | : MS Scan (35 - 450 m/z)                     |
|                      |                                              |
|                      |                                              |

| Į | Handy  | γTI | <u>D</u> ( | <u>Con</u> | ditic | <u>)ns</u> | 5   |   |
|---|--------|-----|------------|------------|-------|------------|-----|---|
|   | Jocorl | ЬT  | -<br>-     | <b>n</b> n |       |            | 200 | 0 |

Desorb Temp.: 300 °C (5 min)Pre Desorb Press.: 180 kPa

#### **VOCs After Combustion of Each Fuel**

We compared the volatile components of burned lumber with gasoline, kerosene, and light oil applied to the lumber. There was a difference in the peak pattern for each fuel.



Chromatogram of volatile components after combustion of each fuel

0.3-0.2-0.1-0.0-

#### Comparison of fuel and volatile components after combustion

It was confirmed that the components of the fuel used for each combustion were detected in the chromatogram after combustion.





### **Related products**

# <Portable Thermal Desorber HandyTD TD265>

Using the GC inlet, the sample is introduced by heat desorption of the volatile components concentrated in the adsorbent.

Due to its compact design that can be carried around, it does not occupy the GC.



| Description   | Qty. | Cat.No     |
|---------------|------|------------|
| HandyTD TD265 | 1    | 2702-30001 |

# <Sampling Bag>

| Description | Capacity | Model | Qty. | Cat.No     |
|-------------|----------|-------|------|------------|
| Tedlar bag  | 10 L     | AAK   | 1 pc | 3008-91110 |

# <Sampling Pump SP209 series>

| Description                              | Qty. | Cat.No     |
|------------------------------------------|------|------------|
| Atmospheric sampling pump SP209-1000Dual | 1 pc | 2702-17583 |

# <Adsorbent for Thermal Desorption AERO TD Tube>

| Description           | Qty. | Cat.No     |
|-----------------------|------|------------|
| AERO TD Tube TENAX TA | 10   | 1003-74102 |

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan 22-1 Nishishinjuku 6-chome Shinjuku-ku, Tokyo 163-1130, Japan

 Phone:
 +81-3-5323-6620

 Fax:
 +81-3-5323-6621

 Email:
 world@gls.co.jp

 Web:
 www.glsciences.com

# www.gisciences.com

#### International Distributors Visit our Website at www.glsciences.com/distributors

GL Sciences Inc. USA 4733 Torrance Blvd. Suite 255 Torrance, CA 90503

USA

Phone: +1-310-265-4424 Fax: +1-310-265-4425 Email: info@glsciencesinc.com Web: www.glsciencesinc.com **GL Sciences B.V.** Dillenburgstraat 7C 5652AM, Eindhoven The Netherlands

Phone: +31-40-254-9531 Email: info@glsciences.eu Web: www.glsciences.eu

#### GL Sciences (Shanghai) Limited

Tower B, Room 2003 Far East International Plaza No.317 Xianxia Road, Changning District Shanghai, China 200051

Phone: +86-21-62782272 Email: contact@glsciences.com.cn Web: www.glsciences.com.cn