

Analysis of Phthalic Acid Effect of Mobile Phase pH on Retention Time

When the pH of mobile phase is changed in the analysis of reverse phase mode on HPLC, the dissociation / non-dissociation state changes depending on the compound, so the retention time changes. Normally, the retention time is longer in non-dissociated state and shorter in dissociated state. This time, using Inertsil WP300 C18, we analyzed the basic compound saccharin and the acidic compound phthalic anhydride standard, and found that the elution order changed under acidic and neutral conditions. We also report that it was estimated that phthalic anhydride in the sample was partially changed to phthalic acid. (R.Takahashi)

Measurement Examples of Saccharin and Phthalic Anhydride



Phosphoric Acid Buffering Action and Phthalic Acid Dissociation



Phosphoric acid has a higher buffering capacity near pH = 2.15, 7.20, 12.35. Even with the same phosphate buffer solution, the retention of acidic compounds is strong under acidic conditions, while the retention of basic compounds is strong under basic conditions.



HPLC condition

Column : Inertsil WP300 C18


5 μm, 150 × 4.6 mm I.D.

Flow rate : 1.0 mL/min Column Temp. : 20 $^{\circ}$ C Detection : UV 230 nm Inj. Vol. : 10 μ L

Gradient Cond. : The solvent is described on the


previous page.

, , , -			
7	Γime(min)	A(vol%)	B(vol%)
	0	10	90
	7.0	10	90
	8.0	95	5
	11.0	95	5



Phthalic anhydride standard product

When the standard phthalic anhydride was measured, two peaks were seen, and it was confirmed that they were phthalic anhydride and phthalic anhydride, respectively.



Comparing phthalic acid and phthalic anhydride, there was a difference in the absorption curves around 270 nm to 280 nm.

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome Shinjuku-ku, Tokyo 163-1130, Japan

Phone: +81-3-5323-6620

Fax: +81-3-5323-6621

Email: world@gls.co.jp

Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255 Torrance, CA 90503

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C 5652AM, Eindhoven The Netherlands

Phone: +31-40-254-9531 Email: info@glsciences.eu Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003 Far East International Plaza No.317 Xianxia Road, Changning District Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

