

# Exosome Purification Columns EVSecond L 70

# Rapid, High purity, Useful Exosome Purification Column



Recent studies have reported significant roles of extracellular vesicle "Exosome" in development and progression of various diseases including cancer metastasis. Therefore exosomes are considered as attractive targets for biomarkers and drug development. However, it remains difficult to isolate high-purity exosomes from biological fluids such as serum.

EVSecond is a size exclusion chromatography open column optimized for effective purification of exosomes. Highly-purified exosomes can be easily collected from serum, plasma, or cell culture supernatant.

EVSecond L70 is a better version of product which the purification can be done at room temperature.



#### **Features**

- Simple gravity-flow handling without ultracentrifugation.
- EVSecond-purified exosomes possess efficient purity for comprehensive miRNA, proteome, and metabolome analysis.
- Exosomes are gently eluted in PBS without structural damage, allowing readministration experiments of collected exosomes to cells or animals.
- EVSecond L70 can used at room temperature and the preparation of column is now much easier.


## **Specification**

	NEW EVSecond L70	EVSecond	
Full length	110 mm	75 mm	
Packing bed length	70 mm	50 mm	
Storing	In a refrigerator at 4 °C	In a refrigerator at 4 °C	
Operation temperature	Room temperature	at 4 °C	
Column preparation	Room temperature Several minutes	at 4 °C Several hours	
Sample volume	50∼1,500 µL	50∼700 µL	

#### [ Operation flow ]


Gravity-flow is applied to each step.

\* Exosome-containing fractions can be identified by western blotting or ELISA experiments detecting tetraspanins (CD9, CD63, CD81, etc.)




# [ Advantages over traditional procedures ]


We compared some purified methods. We analyzed purified exosome using each method in LC/MS, and showed the subcellular localizations of identified proteins. EVSecond provided much higher purification efficiency with less free secreted proteins compared to other existing methods.



### [ Purification of exosomes from human serum ]

Like EVSecond, EVSecond L70 can clearly separate exosomes from serum free proteins.





Isolating exosomes from 200  $\mu$ L of human serum (100  $\mu$ l / fraction)

Red line: CD9-CD9 exosome sandwich ELISA (detecting exosomes)

Blue line: Bradford assay (detecting serum free proteins)

Data provided by Dr. Koji Ueda.

# 【 Ordering Information 】


Description	Quantity	Cat.No.
EVSecond L70 NEW	10 pcs	5010-21395
EVSecond	10 pcs	5010-21390
EvSecond	25 pcs	5010-21392

EVSecond was developed based on the cooperation from Dr. Koji Ueda from Japanese foundation for cancer research.



### [ GL-SPE EXO Fraction Rack ]

This is a fraction rack designed for using EVSecond L70 or EVSecond, which offers smooth column handling and fractionation.



①Column hanger
Maxium 6 columns can be installed.

#### (2)Collection tube racks

Maximum 8 collection tubes (1,5 mL or 2 mL) can be set. The collection tube can be switched by pulling the lever.

#### 3 Base plate

Maximum 6 collection tube racks can be set.

Demension (for EVSecond L70): 300(W)×300(D)×189(H) mm

Demension (for EVSecond): 300(W)×300(D)×150(H) mm

This rack can use both EVSecond L70 and EVSecond by changing column hanger and stand.





#### [ Ordering Information ]

Description	Quantity	Cat.No.
GL-SPE EXO Fraction Rack	1 Set	5010-50450

#### **Contact Information**

**GL Sciences, Inc. USA**4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425

www.glsciencesinc.com email: info@glsciencesinc.com GL Sciences B.V. De Sleutel 9 5652 AS Eindhoven The Netherlands Phone: +31 (0) 40 254 95 31

www.glsciences.eu email: info@glsciences.eu GL Sciences, Inc. Japan 22-1 Nishishinjuku 6-Chome Shinjuku-ku, Tokyo, 163-1130, Japan

Phone: +81-3-5323-6620 Fax: +81-3-5323-6621 www.glsciences.com email: world@gls.co.jp

International Distributors
Visit our Website at
www.glsciences.com/distributors