# MonoSelect C18 for HTS


### Monolithic silica-based column for High-Throughput Screening (HTS)

designed for LC and LC/MS






MonoSelect C18 for HTS is based on monolithic silica technology. The very low flow resistance and the excellent separation performances make it suitable for High-Throughput Screenings (HTS). Due to its metal-free hardware, sharp peaks can be obtained also when analyzing metal chelators. The particular structure of MonoSelect C18 for HTS prevents from clogging and delivers high resistance.



#### **Monolithic Silica**

GL Sciences' monolithic silica consists of precisely controlled  $\mu$ m-sized co-continuous silica skeletons and through-pores. The high external porosity of this structure results in a lower pressure than that of particle-packed columns. Due to the mesopores in the silica skeletons, the surface areas of monolithic silica and particle-packed columns are comparable.




#### **Specifications**

• Bonded Phase: Octadecyl• Max. Operating Temperatur• Max. Operating Pressure: 40 MPa• Carbon Loading• End-capping: Yes	:2-7.5 e: :70°C :7 %
---	----------------------------

#### Why choose MonoSelect C18 for HTS?


MonoSelect C18 for HTS allows for equal or even better separation than C18 3  $\mu$ m silica particle packed columns, and it can be used at both low or high flow rates. The back pressure is only 1/6 of the pressure generated by 1.7  $\mu$ m silica particle packed columns, making it suitable for high throughput analysis.



Contact us for more detailed analytical conditions.


#### Suitable elution time for high-throughput analysis

The faster elution of MonoSelect C18 for HTS compared to other C18 columns and its improved sensitivity make it the best choice for high-throughput analysis.




#### Up to the limits of the fastest equilibration

The equilibration time plays an important role in the result repeatability in gradient analysis. MonoSelect C18 for HTS enables faster equilibration compared to packed columns of the same dimensions, which is desirable for rapid analysis.



#### Superior repeatability during continuous injections

Below data was obtained by performing a continuous injection test on a sample containing 11 compounds and comparing the repeatability of the peak area obtained by each column. The superior repeatability of MonoSelect C18 for HTS is attributed to its fast equilibration.




1 min Analysis Time + 1 min Sample Injection Time (the stabilization is performed during this time) = 2 min

	1	2	3	4	5	6	Ø	8	9	10	1	Ave.
MonoSelect C18 for HTS	2.21%	2.76%	8.89%	4.64%	6.51%	8.39%	3.22%	4.97%	3.17%	4.00%	4.66%	4.86%
Column A	16.0%	13.1%	13.7%	11.0%	7.31%	7.25%	5.80%	6.81%	6.64%	4.88%	7.56%	9.10%
Column B	6.61%	6.45%	9.02%	3.36%	6.49%	7.86%	7.44%	7.62%	6.98%	6.80%	8.96%	7.05%

Contact us for more detailed analytical conditions.


#### Stainless steel and glass double structure to prevent absorption

MonoSelect C18 for HTS structure includes an external stainless hardware and an internal glass tube. The structure does not contain any sintered metal frit, which could cause adsorption of metal chelators. This delivers outstanding sensitivity for such analytes.




#### [Analysis of Fumonisin B1]

Fumonisin B1 is a type of mycotoxin, and the several carboxyl groups of this compound could be a metal chelators. With MonoSelect C18 for HTS, it is possible to obtain a sharp peak even when analyzing this kind of compound.



#### [Analysis of ATP]

ATP contains 3 phosphate groups. Also in this case, good peak shape can be obtained with MonoSelect C18 for HTS.




Contact us for more detailed analytical conditions.

A (%)

B (%)


#### Analysis of proteins

Peaks of macromolecules tend to be distorted with columns packed with particles having small pores. On the other hand, MonoSelect for HTS yields sharp peaks of macromolecules due to the mesopores in the monolithic skeletons.



#### **Superior Durability**

When analyzing compounds contained in physiological samples, proteins must be removed before the analysis in order to avoid column clogging caused by the proteins themselves. However, it is difficult to completely remove all the proteins and this generally leads to the column deterioration over time. The large through-pores of MonoSelect C18 for HTS avoids accumulation of interfering compounds in biological samples. This results in longer repeatability over time compared to particle-packed columns.



Contact us for more detailed analytical conditions.

#### Products

#### MonoSelect C18 for HTS Holder Cartridge Kits

Description	Specifications	Cat.No.		
10mm Holder Cartridge Kit	2.1 mm I.D. x 10 mm cartridge 1PC Holder for 10mm cartridge 1PC	5020-10810		
20mm Holder Cartridge Kit	2.1 mm I.D. x 20 mm cartridge 1PC Holder for 20mm cartridge 1PC	5020-10811		

\* The end- fitting is UP type (Parker style).





10mm Holder Kit

20mm Holder Kit

#### **MonoSelect C18 for HTS Cartridges**


Description	I.D. (mm)	Length (mm)	Cat.No.	Qty.
MonoSelect C18 for HTS Cartridge	2.1	10	5020-10812	1PC
		20	5020-10813	1PC

#### **MonoSelect C18 for HTS Holder**

Description	Length of the Cartridge Applicable (mm)	Cat.No.	Qty.	
MonoSelect C18 for HTS Holder	10mm	5020-10814	1PC	
WONDSEIECT CTO IOL HIS HOIDER	20mm	5020-10815	1PC	

\* The end- fitting is UP type (Parker style).

## 61 Sciences Inc.



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome Shinjuku-ku, Tokyo, 163-1130, Japan Phone: +81-3-5323-6620 Fax: +81-3-5323-6621 Email: world@gls.co.jp Web: www.glsciences.com

#### **International Distributors**

Visit our Website at: https://www.glsciences.com/company/distributor.html

<u>GL Sciences B.V.</u>

Dillenburgstraat 7C

5652 AM Eindhoven

Phone: +31 (0)40 254 95 31

Email: info@glsciences.eu

Web: www.glsciences.eu

The Netherlands

#### GL Sciences (ShangHai) Ltd.

Tower B, Room 2003, Far East International Plaza, NO,317 Xianxia Road, Changning District. Shanghai, China P.C. 200032 Phone: +86 (0)21-6278-2272 Email: <u>contact@glsciences.com.cn</u> Web: www.glsciences.com.cn

#### <u>GL Sciences, Inc. USA</u>

4733 Torrance Blvd. Suite 255 Torrance, CA 90503 Phone: 310-265-4424 Fax: 310-265-4425 Email: <u>info@glsciencesinc.com</u> Web: www.glsciencesinc.com

