

Rapid Solid-Phase Extraction of 136 Pesticides in Water Using Disk Cartridges and Alternative Solvents for GC-MS Analysis

Renato Orsino¹, Rina Nagayasu², Reika Takahara², Manabu Takayanagi², Hiroshi Hayashida²

¹Acore Consumíveis, Brazil, ²GL Sciences Inc., Tokyo, Japan

Introduction

Monitoring of pesticide residues in drinking water is essential for public health and regulatory compliance. Gas chromatography–mass spectrometry (GC-MS) is widely used due to its high sensitivity and selectivity. Prior to analysis, solid-phase extraction (SPE) is typically employed to achieve several hundred-fold concentration. Conventional SPE methods require approximately two hours to complete, with most of the time consumed by the sample loading step. In addition, dichloromethane (DCM)—a commonly used elution solvent—raises safety and environmental concerns due to its high toxicity. This study aimed to develop a faster and safer pretreatment method for GC-MS analysis of pesticides in drinking water by:

- Evaluating whether a disk-type SPE cartridge (EZ Cartridge RP-1) can significantly reduce sample loading time
- Establishing DCM-free elution conditions without sacrificing recovery efficiency

Methods

The analytical method used in this study was based on the official Japanese water quality testing procedures, specifically those described in the guidelines: *"Enactment of Ministerial Ordinance on Water Quality Standards, Partial Revision of Water Supply Law Enforcement Regulations, and Guidelines for Water Quality Management"*—Annex Methods 5 and 5-2. These are authorized procedures for pesticide monitoring in drinking water in Japan.

- SPE Cartridge : EZ Cartridge RP-1 (disk-type, GL Sciences Inc.)
- Target Pesticides : 136 compounds
- Sample Matrix : Ultrapure and tap water spiked at 1 µg/L
- Flow Rates Tested : 15, 50, 100 mL/min
- Elution Solvents:
 - Dichloromethane (DCM)
 - Acetone/Hexane (1:1)
 - Ethyl Acetate




Fig.1 Solid Phase Extraction Procedure

Table 1. Comparison of SPE Cartridge Structure and Sample Flow Rate

	EZ Cartridge	Conventional SPE Cartridge
Structure		
Particle size and diffusion efficiency	10 µm	60 - 70 µm
Sample Volume	500 mL	500 mL
Flow Rate	100 mL/min	10 mL/min
Time	5 min	50 min
Relationship between liner velocity and sample passage area		

Results

Study 1: Effect of Sample Flow Rate on Recovery Rates

Ultrapure water samples were spiked with 136 pesticides at 1 µg/L each, and solid-phase extraction was performed using a disk-type SPE cartridge under three flow rate conditions: 15, 50, and 100 mL/min. The number of compounds recovered with ≥70% efficiency was 131 at 15 mL/min, 132 at 50 mL/min, and 121 at 100 mL/min (Table 3). These results indicate that even at a high flow rate of 100 mL/min, sufficient retention and recovery can be achieved for most compounds, demonstrating that faster flow rates are effective for reducing sample preparation time.




Table 2 GC-MS Conditions

System	GCMS-QP2010 Plus (Shimadzu)
Column	InertCap 5MS/Sil (GL Sciences Inc.) 0.25 mm I.D. × 30 m df = 0.25 µm
Col. Temp.	50 °C (3 min hold) - 10 °C/min - 200 °C - 3 °C/min - 230 °C (5 min hold) - 5 °C/min - 300 °C (2 min hold)
Carrier Gas	He, 100 kPa
Injection	Pulsed Splitless, 1min, 250 °C
Detection	MS SIM
Interface Temp.	280 °C
Sample Size	1.0 µL




Fig.2 Total Ion Current Chromatogram (500µg/L each)

Table 3 Repeatability Linearity, and Recovery (Elution solvent; DCM)

No.	Compounds	Flow Rate		15 mL/min		50 mL/min		100 mL/min		No.	Compounds	Flow Rate		15 mL/min		50 mL/min		100 mL/min	
		Recovery Rate (%)	CV (%)	Recovery Rate (%)	CV (%)	Recovery Rate (%)	CV (%)	Recovery Rate (%)	CV (%)			Recovery Rate (%)	CV (%)	Recovery Rate (%)	CV (%)	Recovery Rate (%)	CV (%)		
1	Dichlorvos	84.4	8.1	82.8	2.9	74.9	6.7	71	Phenothate	85.8	7.3	88.1	7.3	81.5	0.3				
2	Trichlorfon(DEP)	84.4	9.3	82.7	2.1	73	6.7	72	Captan	84.8	4.2	89.2	3.6	80.9	2				
3	Dichlobenil	81.3	6.8	81.8	2.4	73.8	2.6	73	Procymidone	86.4	4.7	92.1	4.3	88.7	1.3				
4	Etridiazole	82.9	10.2	80	4.3	71.1	4.9	74	Dimepiperate	84.4	6.2	89.3	6.9	79.6	0.7				
5	Chloroneb	91.9	5.8	81.3	0.9	83.5	0.1	75	Triflumizole	73.3	10.3	82.2	9.5	66	1.4				
6	Isopropcarb	87.2	5.2	88.9	1.9	83.8	1.8	76	Butamifos Oxon	94.8	11.9	92.6	9.7	75.7	0.7				
7	Molinate	80.6	5.2	82.5	1.2	77.1	3.3	77	Methidathion	87.1	7.8	90.9	5.2	87.6	1.2				
8	Fenobucarb	86.8	4.9	90.4	3	85	1.4	78	Propaphos	80.7	10.1	62.8	7.2	73.7	0.7				
9	Propocur(PHC)	83.4	4.6	92.8	4.6	79.4	1.4	79	Tetrachlorvinphos(CVMP)	83.3	8.6	91.9	8.7	81.4	1				
10	Trifluralin	72.5	4.7	75.8	4.1	58.6	2.5	80	Paclobutrazol	85.9	9.9	91.2	7.6	81.4	0.9				
11	Benfluralin	67.1	6.1	70.8	6.3	53.8	3.6	81	Butachlor	79.1	6	87.2	5.9	77.1	1.9				
12	Cadusafos	79	6.7	83.6	5.4	75.8	4.3	82	alpha-Endosulfan	82.7	4.5	85.8	5.2	76.9	1.2				
13	Pecycuron	88.6	9	95.9	5.7	88.2	0.3	83	9-Bromosanthranicene	-	-	-	-	-	-				
14	Dimethoate	91.1	8.8	96.2	1.7	63.2	2.8	84	Butamifos	89.2	8.5	90.1	8	80.6	1.4				
15	Simazine	87.6	5.2	81.8	5.4	85.3	0.8	85	Napropamide	88.4	6.9	91.7	5.3	89.9	2.8				
16	Atrazine	85.9	4.8	93.9	5.5	86.5	1.5	86	Flutolanil	89.6	9	94.5	5.6	91.3	1				
17	Diazinon Oxon	92.8	8.2	91	7.9	83	2.6	87	(E)-Metominostrubin	86.6	9.2	94.1	5.4	84.2	1.3				
18	Cyanophos(CYAP)	80.5	4.2	89.6	5.6	77.9	0.4	88	Pretlialchor	86.8	7.3	90.3	5.6	86.3	1.4				
19	Propazamide	89.8	7.1	92.9	7.1	87.7	0.4	89	Isoprotiholane	88.3	7.3	91.9	4.5	86.9	2.2				
20	Diazinon	81.1	3.6	87.1	6.2	81.6	2.8	90	Isosathion Oxon	86.5	7.8	77.7	5.5	85.8	0.5				
21	Pyroquilon	84.4	3.9	91.7	3.7	81.7	0.4	91	Uniconazole P	84.8	10.4	91.7	7.4	80.4	0.6				
22	Chlorothalonil(TPN)	74.1	4.1	83.4	5.7	78.7	2.8	92	Thifluzamide	82.3	10.8	92.3	6.7	81.7	1				
23	Anthracene-d10	-	-	-	-	-	-	93	MPP Oxon Sulfoxide	94.5	14.4	150.2	5	84.4	10.4				
24	Ethylthiometon	80.7	4	42.9	3.5	75.9	0.6	94	MPP Oxon Sulfone	91.3	13.8	101.6	6.8	83.4	3.4				
25	Iprobenfos	87.9	6.8	91.1	7.9	79.4	0.6	95	Buprofezin	83.9	6.8	88.7	4.8	78.8	0.4				
26	Tolclofos-methyl Oxon	94.3	8.2	94.4	7.2	84.4	0.4	96	Cyproconazole	83.3	9.2	89.5	7.8	79.9	1.5				
27	Benfuresate	83.9	5.8	93.6	4.7	82.1	2.1	97	Isosathion	86.1	7	89.8	4.9	74.5	1				
28	Dichlofenthion	76	3.7	83.3	5.6	68.7	0.5	98	(Z)-Pyriminobac-methyl	82.2	11.2	92.2	5.8	82.5	0.8				
29	MEP Oxon	97.8	10.6	96.6	8.2	82.6	0.1	99	MPP sulfoxide	93.6	11.3	107.8	3.9	81.4	2.5				
30	Terbutcarb	85.8	6	93.9	6.2	85.9	1.4	100	beta-Endosulfan	85.4	6.2	91.1	4.6	76.5	4.5				
31	Propanil(DCPA)	86.4	5.6	98.4	5.7	86.3	3.7	101	MPP sulfone	92.4	10.8	94.1	6.8	82.5	0.9				
32	Bromobutide	85	6.3	92.3	7.5	84.4	0.8	102	Mespril	91.8	10.9	92.7	3.4	87.8	0.7				
33	Chlorpyrifos-methyl	77.7	4.8	86.8	6.1	75.5	1.4	103	Chlorintrofen(CNP)	90.8	9	88.9	5.4	72.6	4.9				
34	Metribuzin	84.1	7.1	94.8	6.6	80.6	0.8	104	Edifenphos	92	8.9	92.1	4.8	91	0.1				
35	Malaoxon	88.2	15.5	111.6	8.7	89.3	1.8	105	Propiconazole1	90.8	10.4	90.8	6.5	88.4	0.7				
36	Simeconazole	83.1	7.7	94.4	8.7	80.8	1.2	106	Endsulfate	88.2	6.8	91.3	6.9	77.9	4.4				
37	Alachlor	85.4	5	91.8	5.7	85.3	1.1	107	(E)-Pyriminobac-methyl	83.3	10.2	88.6	6.7	84.5	0.9				
38	Tolclofos-methyl	82	3.1	87.7	4.7	79.6	1.7	108	Propiconazole2	86.9	9.6	90.3	6.8	68.5	0.9				
39	Simetryne	80.8	7.4	84.1	9.4	82.2	0.3	109	EPN Oxon	99.1	11.9	101.7	5.8	85.8	3.1				
40	Metaxyl	86.9	7.2	95.5	6.1	85.2	1.7	110	Thenylchlor	92.6	8.6	91.7	4.7	90.4	1.1				
41	Ametryn	79.8	6.9	91.1	8.3	81.7	1.7	111	Tebuconazole	85.3	9.4	93.9	6.7	80.9	0.7				
42	Cimethylin	78.9	3.5	88.4	4.1	80.2	0.4	112	Propargite(BPPS)	59.4	9.3	77.7	7.1	41.3	1.3				
43	Prometryn	77.3	6.8	90.8	7.4	79.1	0.6	113	Pyributicarb	82.9	10.7	77.8	7.9	60.6	0.4				
44	Dithiopyr	73.4	4.7	83.6	7.8	67.3	2	114	Pyridaphenthion	87.5	14.5	89.6	8.7	79.3	0.1				
45	MPP Oxon	90.7	8.5	73.6	6.5	80.1	0.7	115	Acetamidip	84.7	14.4	91.4	11.1	70.6	6.1				
46	Pirimiphos-methyl	76	4.7	86.9	6.8	75.1	0.8	116	Iprodion	86.2	10.1	90.3	5.1	83.5	0.1				
47	Fenitrothion	91.2	8.3	94.4	8.4	83.6	1.9	117	Chrysene-d12	-	-	-	-	-	-				
48	Bromacil	86.2	11.2	100.5	9.6	81.4	0.6	118	EPN	91.2	9.9	90.2	6.5	75	1.4				
49	(E)-Dimethylvinphos	84	6.7	98.7	9.9	83.9	1.4	119	Piperophos	82.6	12.4	84.6	10.3	72.5	0.4				
50	Esprocarb	82.7	4	89.6	4.4	82.1	2	120	Cumyluron	81.5	11.5	88.4	7.5	69.3	3.4				
51	Malathion	86.8	6.1	95.2	6.4	87.5	1.6	121	Indanofan	83.1	11.1	54.2	8.4	75	1.7				
52	Chlorpyrifos Oxon	93.1	11.6	101.1	10.2	86.7	1.2	122	Anilofos	85.8	13.2	88.2	9.3	76.7	0.2				
53	Quinoclamine(ACN)	67.9	5.5	96.5	5.5	76.7	2.1	123	Oryastrobin	83.5	12.9	89.7	8.4	74.5	1.6				
54	Metolachlor	78.6	4.4	91.9	6.4	81.3	2.1	124	Bifenox	94	16.7	91.1	8.5	68.7	3.4				
55	Chlorpyrifos	77.8	4.1	85.4	6.9	71	2	125	Furametryr	82.4	11.1	89.6	6.2	81.4	0.4				
56	Thiobencarb	82.6	3.2	87.6	4	85	1.8	126	Iprodion-t	93.2	9.3	89.6	5.6	74.4	1.3				
57	(Z)-Dimethylvinphos	82.7	8.2	91.2	8.3	82.7	1.5	127	Phosalone	80.8	10.7	87.5	7.7	74	0.3				
58	Cyanazine	83.2	7.6	91.4	6.5	82.3	1.4	128	Pyriproxyfen	77.9	7.6	82.9	4.4	64.2	0.9				
59	Fenthion	83.4	3.8	75.9	5.3	82.4	2.2	129	Mefenacat	89.3	10	91.4	6.4	81.9	0				
60	Chlorthal-dimethyl(TCTP)	79.2	2.8	85.2	5.1	77	3.2	130	Cyhalofop-butyl	55.5	9.3	70.4	4.4	44.7	0.4				
61	Isofenphos Oxon	95.3	13	98.6	10.9	81.5	0.1	131	CPA-amino	92	11.1	92.4	2.4	85.5	3.3				
62	Tetraconazole	80.3	9.3	91.5	8.7	80.9	1.6	132	Pyradofos	174.9	50.8	83.8	48.1	128.1	2				
63	Fthalide	86.5	4.6	90.7	3	84.9	4.6	133	Etofenbazine	93.5	14.9	88.6	5.6	82.3	1.8				
64	Fosfthiazate	86.2	9.7	105.4	6.9	85.5	0.2	134	Castenrolle	97.4	12	91.8	5.6	91	0.5				
65	Thiamethoxam	70.1	15.2	88.7	2.7	39.3	2.5	135	Boscalid	95.6	13	87.4	3.5	86.7	1				
66	Pendimethalin	79.6	5.7	81.7	7.4	63.5	2	136	Ethofenprox	94.5	10.8	92.9	4	35.4	0.6				
67	Cyprodinil	79.8	5.7	81.7	6.2	79.2	0.3	137	Thiacloprid	51.3	10.8	49.6	2.1	79.5	7				
68	Isofenphos	82.7	8.5	87.5	8.3	84.7	0.5	138	Difenoconazole	89.6	15.2	91.3	4.7	70.8	4.2				
69	Isofenphos	82.9	6.1	87.8	4.8	82.6	1	139	Pyraoxifen	91.6	9.9	91.3	4.7	72.3	1.3				
70	Methyldymron	86	6.4	91	4.5	82.6	1.1												