

Concentration Analysis of Antiperspirant Flavor Components - Using HandyTD Portable Thermal Desorber and MonoTrap Simplified Enrichment Tools

The flavor components of liquid antiperspirants were screened and analyzed using HandyTD and MonoTrap RGPS TD a simplified enrichment tool. HandyTD is a compact, portable, and easy-to-operate thermal desorption device. After sample collection, set the MonoTrap in a glass liner and introduce the heat-desorbed components into a GC inlet. Terpenes and other substances that confer refreshment were detected.

Preliminary processing procedure

Liquid antiperspirant

Place 500 μ L of sample in a 40 mL vialCollection (HS)
One MonoTrap RGPS TD

15 min at room temperature

TD/GC/MS




GC Conditions

System	: Thermal Desorption-GC/MS (HandyTD TD265)
Column	: InertCap Pure-WAX 0.32 mm I.D. x 60 m, df = 0.50 μ m
Col. Cat. No.	: 1010-68264
Col.Temp.	: 40 °C(5 min) - 6 °C/min - 250 °C
Carrier Gas	: He, 2.0 mL/min
GC Inlet	: 250 °C, Split 9:1
Detection	: MS Scan (m/z 45-450)

HandyTD Conditions

Desorb Temp.	: 40 °C - 45 °C/sec - 250 °C (1.5 min)
Desorb Press.	: 114 kPa



1. Ethanol
2. Pinene
3. Pinene
4. Pinene

5. Myrcene
6. Terpinolene
7. Limonene
8. Terpinyl methyl ether

9. Cymene
10. Linalool
11. Isopulegol
12. Levomenthol

* Standard samples are not used for qualitative analysis. Results from a library search.

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsciences.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO,317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

International Distributors

Visit our Website at www.glsciences.com/distributors