GC Technical Note

GT141 GL Sciences Inc.

Analysis of Ethanol in Food

Presented here is an analysis of ethanol in food using FID mounted gas chromatograph dedicated to packed column. For the analysis of sample containing non-volatile or refractory components, selecting an insert injection method and by keeping the residual contaminants within an injection port liner, the impact on analytical column can be reduced. In the analysis of soy sauce or soybean paste used in this example, the sample can be quantitated only by dilution.

Analysis Example

<Standard Sample>



1. Methanol 2. Ethanol 3. 2-Propanol 4. 1-Propanol (each 0.5 % in H2O, v/v)

Conditions Column Col.Temp. Carrier Gas Detector Injection Sample Size

:Gaskuropack 54 60/80 Glass 1/4" O.D. X 2 m X 3.0 mml.D. :150 °C : N_2 30 mL/min : FID 200 °C : 200 °C : 1 μ L

<Actual Sample>



[Soybean paste] ••• Diluted to 25mg/mL by purified water



Ethanol contained in soy sauce and soybean paste was analyzed. Each sample was diluted by an arbitrary amount of purified water and centrifuged (5,000xg, 1 min.). For 1mL of supernatant, 10 μL of 50% internal standard material was added. Internal reference method was applied as an analysis method and 2-propanol was used as internal standard material.

And an insert injection method was adopted as injection method.



Insert Injection Port



When glass column is used, use an insert injection port so that non-volatile or refractory component may not be introduced directly into the column. Photo on the left-hand side shows an injection port liner right after injection of soy sauce and soybean paste. It is known that non-volatile component contained in the sample remained at the injection port liner after usage. By using an insert injection port, particularly useful results can be obtained, enabling the quantitation only by dilution of sample as well as the extension of column life time.





When an injection port liner is contaminated, replace with a new one.

GL Sciences B.V.

5652 AS Eindhoven

Phone: +31 (0)40 254 95 31

Email: info@glsciences.eu

Web: www.glsciences.eu

The Netherlands

De Sleutel 9

Related Products

Column

- Packing material Gaskuropack 54 60/80 20 g Cat.No.1002-45406
- Empty column Glass 1/4" O.D. x 2 m x 3.0 mml.D. Cat.No.3003-64120
- * Packed column also can be manufactured.

- Injection port liner
 - Glass liner(for GC3220, GC 3210, andGC-4000 Plus) Cat.No. 2701-22753
- ●Wool
 - Quartz wool Fine(1 $^{\sim}$ 6 $\mu m)$ Silane finish Cat.No. 3001-12404

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan 22-1 Nishishinjuku 6-Chome Shinjuku-ku, Tokyo, 163-1130, Japan Phone: +81-3-5323-6620 Fax: +81-3-5323-6621 Email: world@gls.co.jp Web: www.glsciences.com

International Distributors Visit our Website at:

https://www.glsciences.com/company/distributor.html

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003, Far East International Plaza, NO,317 Xianxia Road, Changning District. Shanghai, China P.C. 200032 Phone: +86 (0)21-6278-2272 Email: <u>contact@glsciences.com.cn</u> Web: www.glsciences.com.cn GL Sciences, Inc. USA 4733 Torrance Blvd. Suite 255 Torrance, CA 90503 Phone: 310-265-4424 Fax: 310-265-4425 Email: info@glsciencesinc.com Web: www.glsciencesinc.com

