


## Analysis of Ascorbic acid [Vitamin C] (2)



<b>System</b>	: GL-7400 system	<b>Solutes</b>	: 1. Erythorbic acid (5 mg/L)
<b>Column</b>	: Inertsil NH2 (5 μm, 250 x 4.6 mm I.D.)		2. Ascorbic acid (5 mg/L)
<b>Column Cat. No.</b>	:		
<b>Eluent</b>	: A) CH <sub>3</sub> CN B) CH <sub>3</sub> OH C) 10 mM NaH <sub>2</sub> PO <sub>4</sub> D) 0.03 % Homocysteine A/B/C/D = 600/30/100/30 ,v/v/v/v		
<b>Flow rate</b>	: 1.0 mL/min		
<b>Col. Temp.</b>	: 40 °C		
<b>Detection</b>	: PDA 270 nm (GL-7452A PDA Detector)		
<b>Injection Vol.</b>	: 5 μL		
<b>Sample</b>	:		