


## Analysis of Vitamin E



### Conditions

**System** : GL-7400 HPLC system  
**Column** : InertSustain NH2  
 (5 $\mu$ m, 250 x 4.6 mm I.D.)  
**Column Cat. No.** : 5020-16628  
**Eluent** : A) *n*-Hexane  
 B) Ethyl acetate  
 A/B = 70/30, v/v  
**Flow Rate** : 1.0 mL/min  
**Col. Temp.** : 30 °C  
**Detection** : UV 290 nm (GL-7450 UV Detector)  
**Injection Vol.** : 10  $\mu$ L

### Analyte:

1. $\alpha$ -Tocopherol	25 mg/L
2. $\beta$ -Tocopherol	25 mg/L
3. $\gamma$ -Tocopherol	25 mg/L
4. $\delta$ -Tocopherol	25 mg/L