


Analysis of Haloperidol

(Under the Condition of the Japanese Pharmacopoeia 17th edition)



Conditions

System	: GL7700 HPLC system
Column	: Inertsil ODS-4 (5 μ m, 150 x 4.6 mm I.D.)
Column Cat. No.	: 5020-03945
Eluent	: Buffer*
Flow Rate	: 1.0 mL/min
Col. Temp.	: 40 °C
Detection	: UV 220 nm (UV7750 UV Detector)
Injection Vol.	: 10 μ L
Sample	: Standard

Analyte

1. Haloperidol

The peak area ratio of 1-a to 1-b : (15 \leq) 20.0 (\leq 25)

Theoretical plates of 1-b	: 9,176 (\geq 4,000)
Tailing factor of 1-b	: 1.06 (\leq 2.0)
RSD of the	
1-b area (%) (n=6)	: 0.55 (\leq 2.0)

*Dissolve 2.95 g of trisodium citrate dihydrate in 900 mL of water.

Adjust to pH 3.5 with diluted hydrochloric acid.

Add water to make 1,000 mL.(Solution A)

Add 700 mL of CH₃OH to 300 mL of solution A.(Solution B)

Add 1.0 g of sodium dodecyl sulfate to Solution B.