


Analysis of Mirtazapine

(Under the Condition of draft for USP, Mirtazapine Compounded Oral Suspension, Veterinary)



Conditions

System	: GL7700 HPLC system	Analyte:	
Column	: Inertsil ODS-3 (5 μ m, 250 x 4.6 mm I.D.)	1. Mirtazapine	500 mg/L
Column Cat. No.	: 5020-01732		
Eluent	: Buffer*		
Flow rate	: 1.5 mL/min		
Col. Temp.	: 60 °C	Symmetry factor (1)	: 1.06 (< 2.0)
Detection	: UV 295 nm (UV7750 UV Detector)	RSD of the peak area of 1 (%) (n=6)	: 0.18 (< 2.0)
Injection Vol.	: 10 μ L		
Sample	: Standard		

*Add 23.5 mL of 25% tetramethylammonium hydroxide solution to 626.5 mL of water.

Adjust with phosphoric acid to a pH of 7.4. Add 150 mL of acetonitrile, 125 mL of methanol, and 75 mL of tetrahydrofuran.