


# Analysis of Etoposide phosphate

(Under the Condition of draft for USP)



## Conditions

**System** : GL7700 HPLC system  
**Column** : InertSustain Phenyl  
           (5  $\mu$  m, 150 x 4.6 mm I.D.)

**Column Cat. No. :** 5020-16327

**Eluent** : A) CH<sub>3</sub>CN  
           B) SolutionA\*  
           A/B = 15/85, v/v\*\*

**Flow Rate** : 1.5 mL/min

**Col. Temp.** : 40 °C

**Detection** : UV 230 nm (PD7752 PDA Detector)

**Injection Vol.** : 25  $\mu$  L

**Sample** : Standard

\* 2.3g/L of ammonium phosphate, monobasic in water.

\*\* Acetonitrile and SolutionA (15:85). Adjust with phosphoric acid to a pH 2.5.

## Analyte

1. Etoposide phosphate  
 2. Etoposide

RSD of the peak area of 1-a (%)(n=6)

: 2.77 ( $\leq$  10)

RSD of the peak area of 2 (%)(n=6)

: 2.85 ( $\leq$  10)

S/N for 1-b

: 41.3 ( $\geq$  10)