


# Analysis of Sucrose

(Under the condition of the PDG (Pharmacopeial Discussion Group) Stage 2 Drafts)



## Conditions

System	: Chromaster
Column	: InertSphere Sugar-2 (9 $\mu$ m, 300 x 7.8 mmI.D.)
Column Cat. No.	: 5020-11000
Eluent	: H <sub>2</sub> O
Flow Rate	: 0.5 mL/min
Col. Temp.	: 80 °C
Detection	: RI (5450 RI Detector)
Injection Vol.	: 10 $\mu$ L
Sample	: Standard

## Analyte:

1. Raffinose : Relative retention time (1 to 2-a) : 0.92 (about 0.9)
2. Sucrose
3. Glucose : Relative retention time (3 to 2-a) : 1.19 (about 1.2)
4. Fructose : Relative retention time (4 to 2-a) : 1.48 (about 1.6)

## System suitability solution

- Peak-to-valley ratio (1, 2-a) : 31.07 ( $\geq 2.5$ )  
Resolution (2-a, 3) : 3.86 ( $\geq 1.5$ )

## Impurity Standard solution

- RSD of the Peak area of 1 (%) (n=6) : 1.16 ( $\leq 5.0$ )  
RSD of the Peak area of 2-b (%) (n=6) : 1.76 ( $\leq 5.0$ )  
RSD of the Peak area of 3 (%) (n=6) : 0.61 ( $\leq 5.0$ )  
RSD of the Peak area of 4 (%) (n=6) : 1.19 ( $\leq 5.0$ )