


Analysis of Acyclovir

(Under the Condition of the Japanese Pharmacopoeia 18th, Acyclovir Injection)



Conditions

System	:	Chromaster PLUS HPLC system (HITACHI)	
Column	:	InertSustain C18 (GL Sciences Inc.) (5 μ m, 250 x 4.6 mm I.D.)	1. Nicotinic acid 15 μ g/mL 2. Acyclovir 10 μ g/mL
Column Cat. No.	:	5020-07346	
Eluent	:	CH ₃ OH/Solution*=5/95 ,v/v	Resolution : 11.9 (≥ 3)
Flow Rate	:	0.78 mL/min	
Col. Temp.	:	25 °C	RSD of the peak area ratio of
Detection	:	UV 254 nm	2 to 1 (%)(n=6) : 0.10 (≤ 1.0)
Injection Vol.	:	20 μ L	
Sample	:	Standard	

*: To 1.45 g of phosphoric acid and 25 mL of dilute acetic acid add water to make 900 mL. Adjust this solution to pH 2.5 with 1 mol/L sodium hydroxide, and add water to make 1000 mL.