



Analysis of Linagliptin

(Under the Condition of draft for USP PF52(3), Linagliptin, ORGANIC IMPURITIES, PROCEDURE 2)



Sensitivity solution



Standard solution

Conditions

System : Chromaster PLUS (HITACHI)
Column : InertSustain C18 (GL Sciences Inc.)
 (5 μ m, 250 x 4.6 mm I.D.)
Column Cat. No. : 5020-07346
Eluent : A) CH₃OH/CH₃CN/H₂O = 70/10/20, v/v/v
 B) CH₃OH/Buffer* = 10/90, v/v

Analyte:

1-a. Linagliptin 0.25 μ g/mL
 1-b. Linagliptin 0.50 μ g/mL

Signal-to-noise ratio (1-a) : 15.1 ($\geq$ 10)

RSD of
 the peak area of 1-b (%) (n=6) : 0.56 ($\leq$ 5.0)

Time (min)	A (vol %)	B (vol %)
0.0	30	70
12.0	50	50
45.0	85	15
55.0	85	15
60.0	30	70
70.0	30	70

Flow Rate : 1.0 mL/min
Col. Temp. : 45 $^{\circ}$ C
Detection : UV 225 nm (Chromaster 5410AD UV)
Injection Vol. : 20 μ L
Sample : Standard in CH₃OH/ H₂O = 40/60, v/v

*: Dissolve 2.04 g of potassium phosphate, monobasic in 1000 mL of water. Adjust with 10% H₃PO₄ in H₂O to a pH of 3.2.