

Simultaneous Analysis of Cyanoguanidine and Melamine in Food by LC/MS/MS


It became a serious problem in China and Taiwan in 2013 that cyanoguanidine, or dicyandiamide, was detected in powdered milk made in New Zealand. Since cyanoguanidine is contained in some chemical fertilizer as nitrification inhibitor, it is considered as a possible cause of the problem that such fertilizer was used in ranch.

In 2008, melamine was detected in milk and dairy products, which also became a big issue. HILIC



(hydrophilic interaction chromatography) mode was used for the determination of melamine content because melamine is highly hydrophilic. Chemical structure of cyanoguanidine is similar to that of melamine. In this note, cyanoguanidine was analyzed with HILIC mode. Cyanoguanidine and melamine could be determined simultaneously using Inertsil HILIC as a separation column.

(K.Kanno)

Chromatogram Obtained from Standard Sample



Chemical Structure



Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.

HPLC Condition


System : LC800
Column : Inertsil HILIC
(5 µm, 150 x 3.0 mm I.D.)
Eluent : A) CH₃CN
B) 10 mM Ammonium acetate
A/B = 90/10 — 0.5 min — 90/10 — 5.5 min —
50/50 (Equilibration for 5 min), v/v
Flow Rate : 0.5 mL/min
Col. Temp. : 40 °C
Detection : MRM (SRM)
Inj. Vol. : 5.0 µL

MS/MS (MRM) Condition

System : 4000 QTRAP (ABSciex)
Detection : Compound

	Q1	Q3	DP	CE	CXP	EP
Cyanoguanidine	85.1	68.1	41	41	6	10
Melamine	126.9	68.1	41	41	10	10

Ion Source : ESI (Posi) CUR CAD IS TEM GS1 GS2 ihe
10 8 5500 700 60 30 on



Calibration Curve

Chromatograms Obtained from Powdered Milk Sample

Example of
Sample Pretreatment


Sample

1.0 g


Extraction

- 20 mL of water/acetonitrile = 50/50 (v/v)
- Sonication extraction for 30 min
- Filtration (0.45 µm membrane filter)

HPLC



1. Cyanoguanidine
2. Melamine



Cyanoguanidine and melamine was added to commercially available powdered milk. Each concentration was adjusted to 10 µg/L (= 0.2 mg/kg powdered milk).

< HPLC Column used in This Note >

Inertsil HILIC 5 µm, 150 x 3.0 mm I.D. (Cat.No. 5020-07715)

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO,317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

International Distributors

Visit our Website at:

<https://www.glsciences.com/company/distributor.html>