

# HPLC Analysis of Vitamin B1, B2


Takadiastase (*Aspergillus oryzae* carboxyl proteinase) has been used for hydrolyses of thiamin phosphate ester and riboflavin phosphate ester to analyse thiamin and riboflavin, respectively by HPLC.

However, it is difficult to get commercial Takadiastase stable nowadays.

So in this technical note, we introduce a new sample preparation method for thiamin and riboflavin analyses using acid phosphatase hydrolysis instead of Takadiastase hydrolysis.

## Standard Solution (Vitamin B1)

Separated thiamine by ODS column was derivatized with potassium ferricyanide under alkaline condition. Fluorescent compound thiochrome was produced and detected by fluorescence detector.



### HPLC Conditions (Vitamin B1)

**Column:** Inertsil ODS-3 5 mm, 150 × 4.6 mm I.D.

**Eluent :** A) CH<sub>3</sub>OH

B) 10 mM NaH<sub>2</sub>PO<sub>4</sub>, 150 mM NaClO<sub>4</sub> (pH 2.2) A/B = 1/9, v/v (Mixed by a gradient mixer)

**Flow Rate:** 1.0 mL/min

**Col. Temp.:** 40 °C

**Detection :** FL Ex. 375 nm, Em. 440 nm


**Inj. Vol. :** 20 µL

**Reaction Reagent:** 0.05 w/v % C<sub>6</sub>FeK<sub>3</sub>N<sub>6</sub>, 15 w/v % NaOH

**Flow Rate of Reaction Reagent:** 0.4 mL/min

**Reaction Coil:** 5000 × 0.33 mm I.D. (PTFE tube)

## Standard Solution (Vitamin B2)



### HPLC Conditions (Vitamin B2)

**Column :** Inertsil ODS-3 5 µm, 150 × 4.6 mm I.D.

**Eluent :** A) CH<sub>3</sub>OH

B) 40 mL of 4 M CH<sub>3</sub>COONa, 20 mL of 50 % CH<sub>3</sub>COOH, 1940 mL of H<sub>2</sub>O (pH 4.5)


A/B = 35/65, v/v (Mixed by a gradient mixer)

**Flow Rate :** 1.0 mL/min


**Col. Temp. :** 40 °C

**Detection :** FL Ex. 445 nm, Em. 530 nm

**Inj. Vol. :** 20 µL




Calibration Curve

**Sample Preparation Method Example****Solid Phase Extraction**

\*1 Please do these steps as necessary.


**An Effect of Solid Phase Extraction (a case of commercial instant noodle)****1. Thiamine**

No Sample preparation



Clean Up


With Sample preparation




## Vitamin B1 in foods

1. Thiamine
2. Hydroxyethylthiamine (HET)

### An extract of beef




### A commercial instant noodle




## Vitamin B2 in foods

### An extract of beef

3. Riboflavin



### A commercial instant noodle




Post column derivatization HPLC system was used for Vitamin B1 analysis. On the other hand, Vitamin B2 is natural fluorescence compound and was analyzed by normal HPLC system. In this time, HPLC system in which switching valve was mounted like below.

## Schematic Diagram

Switching Valve for Change Flow Channels

HPLC Column



### HPLC Column

Inertsil ODS-3 5  $\mu$ m, 150  $\times$  4.6 mm I.D.  
Cat. No. 5020-01731

### SPE Column

InertSep SCX-2 500 mg/6 mL  
Cat. No. 5010-61724

### Syringe Filter

GL Chromato Disk Sample Filters 25A  
PORE SIZE 0.45  $\mu$ m  
Cat. No. 5040-28512

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

### GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome  
Shinjuku-ku, Tokyo,  
163-1130, Japan  
Phone: +81-3-5323-6620  
Fax: +81-3-5323-6621  
Email: [world@glsc.co.jp](mailto:world@glsc.co.jp)  
Web: [www.glsciences.com](http://www.glsciences.com)

### GL Sciences B.V.

De Sleutel 9  
5652 AS Eindhoven  
The Netherlands  
Phone: +31 (0)40 254 95 31  
Email: [info@glsciences.eu](mailto:info@glsciences.eu)  
Web: [www.glsciences.eu](http://www.glsciences.eu)

### GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,  
Far East International Plaza,  
NO,317 Xianxia Road,  
Changning District.  
Shanghai, China P.C. 200032  
Phone: +86 (0)21-6278-2272  
Email: [contact@glsciences.com.cn](mailto:contact@glsciences.com.cn)  
Web: [www.glsciences.com.cn](http://www.glsciences.com.cn)

### GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255  
Torrance, CA 90503  
Phone: 310-265-4424  
Fax: 310-265-4425  
Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)  
Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

### International Distributors

Visit our Website at:

<https://www.glsciences.com/company/distributor.html>