



# Analysis of Steviol Glycoside

## -- Japanese Pharmaceutical Excipients

The Ministry of Health, Labour and Welfare issued a notification on March 29, 2018, for revise the Standards for Japanese Pharmaceutical Excipients 2018 (Notification No. 0329-1 of the Crude Drug Administration). In this Pharmaceutical Excipients Specification 2018, the specifications for the five newly reviewed products were listed, partially revised for the 25 already listed products, and the four ingredients of the eight products were combined, resulting in 489 listed products. This paper presents the results of HPLC analyses of "purified Stevia extract" among the five newly listed products.

(M. Mano)

### Purified Stevia Extract Purity Test (4) Steviol



#### HPLC Conditions

System : GL7700 HPLC system  
 Column : Inertsil ODS-HL  
           (5  $\mu$ m, 250 x 4.6 mm I.D.)  
 Eluent : A) CH<sub>3</sub>CN  
           B) 0.1 % H<sub>3</sub>PO<sub>4</sub> in H<sub>2</sub>O  
           A/B = 55/45, v/v  
 Col. Temp. : 40 °C  
 Detector : UV 213 nm  
 Injection Vol. : 20  $\mu$ L  
 Flow Rate : 1.25 mL/min

#### 【System Suitability Test】


◆ System Stability

(1, 2-a) : 17.5 ( $\geq 6$ )

◆ System Reproducibility

Steviol Peak Area (2-b)

RSD (%) (n=6) : 0.10 ( $\leq 1.5$ )

**Purified Stevia Extract Assay****HPLC Conditions**

System : GL7700 HPLC system  
 Column : Inertsil NH2  
           (5  $\mu$ m, 250 x 4.6 mm I.D.)  
 Eluent : A) CH<sub>3</sub>CN  
           B) H<sub>2</sub>O  
           A/B = 80/20, v/v  
 Col. Temp. : 35 °C  
 Detector : UV 213 nm  
 Injection Vol. : 20  $\mu$ L  
 Flow Rate : 1.40 mL/min

**【System Suitability Test】**


◆ System Stability

**Separation (3-a, 4) : 3.73 ( $\geq 1.5$ )**


◆ System Reproducibility

**Steviol Peak Area (3-b)****RSD (%) (n=6) : 0.08 ( $\leq 1.5$ )**

## Chemical Structure



4-*tert*-Butylbenzoic acid



Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.

Compounds	R <sub>1</sub>	R <sub>2</sub>
Steviol	H	H
Rubusoside	β-Glucose	β-Glucose
Dulcoside A	β-Glucose	β-Glucose - α-Rhamnose (2→1)
Stevioside	β-Glucose - β-Glucose (2→1)	β-Glucose - β-Glucose (2→1)
Rebaudioside C	β-Glucose	β-Glucose - α-Rhamnose (2→1)   β-Glucose (3→1)
Rebaudioside A	β-Glucose	β-Glucose - β-Glucose (2→1)   β-Glucose (3→1)

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

### GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome  
Shinjuku-ku, Tokyo,  
163-1130, Japan  
Phone: +81-3-5323-6620  
Fax: +81-3-5323-6621  
Email: [world@glsciences.co.jp](mailto:world@glsciences.co.jp)  
Web: [www.glsciences.com](http://www.glsciences.com)

### GL Sciences B.V.

De Sleutel 9  
5652 AS Eindhoven  
The Netherlands  
Phone: +31 (0)40 254 95 31  
Email: [info@glsciences.eu](mailto:info@glsciences.eu)  
Web: [www.glsciences.eu](http://www.glsciences.eu)

### GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,  
Far East International Plaza,  
NO,317 Xianxia Road,  
Changning District.  
Shanghai, China P.C. 200032  
Phone: +86 (0)21-6278-2272  
Email: [contact@glsciences.com.cn](mailto:contact@glsciences.com.cn)  
Web: [www.glsciences.com.cn](http://www.glsciences.com.cn)

### GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255  
Torrance, CA 90503  
Phone: 310-265-4424  
Fax: 310-265-4425  
Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)  
Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

### International Distributors

Visit our Website at:

<https://www.glsciences.com/company/distributor.html>