GL Sciences Inc.

Anaysis of Antioxidant capacity by pre-column HPLC using DPPH radical

This time, we will introduce an example of measuring the antioxidant capacity (DPPH radical scavenging activity) of food using 1,1-diphenyl-2-picrylhydrazyl free radical (DPPH), which is known as a stable radical. DPPH is known as a purple stable radical and becomes colorless when reduced. This time, by adding tea to DPPH, we confirmed whether changes would occur due to antioxidants such as catechins.

(K. Suzuki)



Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.



HPLC conditions

Column	: InertSustainSwift C8
	(5 μm, 150 × 4.6 mm I.D.)
Cat. No.	: 5020-88331
Eluent	: A) CH₃CN
	B) H ₂ O
	A/B = 55/45, v/v
Flow rate	: 1.0 mL/min
Col.Temp.	: 40 °C
Detection	: VIS 540 nm (DAD)
Inj.Vol.	: 10 μL



DPPH free radical calibration curve

Analysis example

The sample was added so that the DPPH was 100 mg/L, and the measurement was performed.



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan 22-1 Nishishinjuku 6-chome Shinjuku-ku, Tokyo 163-1130, Japan

 Phone:
 +81-3-5323-6620

 Fax:
 +81-3-5323-6621

 Email:
 world@gls.co.jp

 Web:
 www.glsciences.com

6L Sciences

International Distributors Visit our Website at www.glsciences.com/distributors

GL Sciences Inc. USA 4733 Torrance Blvd. Suite 255 Torrance, CA 90503 USA

 Phone:
 +1-310-265-4424

 Fax:
 +1-310-265-4425

 Email:
 info@glsciencesinc.com

 Web:
 www.glsciencesinc.com

<u>GL Sciences B.V.</u> Dillenburgstraat 7C 5652AM, Eindhoven The Netherlands

Phone: +31-40-254-9531 Email: info@glsciences.eu Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003 Far East International Plaza No.317 Xianxia Road, Changning District Shanghai, China 200051

Phone: +86-21-62782272 Email: contact@glsciences.com.cn Web: www.glsciences.com.cn