


Aflatoxins are mycotoxins produced by molds. HPLC methods are used to analyze them and to reduce the impact of contaminants, pretreatment is required for sample cleanup.

InertSep VRA series for food analysis is a multifunctional Mini-Column (Clean-Up Cartridge) useful for cleanup of mycotoxins and other compounds. VRA-1 is a useful phase for common cereals, and VRA-3 reverse-phase-enriched phase is suited to oil-rich beans.

1. Flow diagram of solid phase pretreatment



Reference: Ministry of Health, Labour and Welfare Food Safety Notification No. 0816-2

2. Products for solid-phase extraction

[InertSep VRAs Series]


This is a multifunctional (multi-function) solid-phase minicolumn dedicated to clean-up processing of organic compounds from high-matrix samples.

Syringe barrel type cartridge

Product name	Detailed	Column Size	Qty.	Cat.No.
InertSep VRA-1	Mixed mode of reversed phase + ion exchange layer	6mL	30 bottles	5010-68140
InertSep VRA-2	Economy version of VRA-1	6mL	30 bottles	5010-68141
InertSep VRA-3	Recommendation for this study Reversed-phase-enriched form of VRA-1	6mL	30 bottles	5010-68142

[Gravity Flow Manifold]

This is a Free fall manifold that can use multifunctional mini column for aflatoxin pretreatment and 4 mL vial as a receptacle.



Product name	Qty.	Cat.No.
GL-SPE Gravity Flow Manifold	Formula 1	5010-50430
GL-SPE 4 mL vial tray for Gravity Flow Manifold	One unit	5010-50432
Luer stop valve (made of PTFE)	12	5010-60010

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO,317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

International Distributors

Visit our Website at www.glsciences.com/distributors