
14. Purification of essential fatty acid in Rapeseed oil



Column	: Inertsil PREP-ODS
	10µm
	250 × 20 mm I.D.
Column Cat. No.	: 5020-31123
Eluent	: CH_3OH
Flow rate	: 1.0 mL/min.
Col. Temp.	: 40°C
Detection	: RI
Injection Vol.	: 3 mL
Sample	: Methyl Esters of Rape
	Seed Oil (0.1 g/mL)

Solutes : 1. γ - Linolenic acid
2. Linolenic acid
3. Oleic acid
4. Stearic acid
5. Arachidic acid