
Drug



Column : Inertsil ODS-3V
150 × 4.6 mm I.D.

Column Cat. No. :

Eluent : CH3CN/0.1 % H_3PO_4
+ 5 mM IPCC-05 = 15/85

Flow rate : 1.0 mL/min.

Col. Temp. : 40ºC

Detection : UV 210 nm

Injection Vol. : 20 µL

Sample : Water-soluble Vitamins (Drugs)

Solutes : 1. Pyridoxine (0.10 mg/mL)
2. Riboflavin (0.10 mg/mL)
3. Fursultiamine (1.00 mg/mL)