
Vitamins



Column : Inertsil CN-3
5µm
250 × 4.6 mm I.D.

Column Cat. No.:

Eluent : CH_3CN/H_2O = 60/40

Flow rate : 1.0 mL/min.

Col. Temp. : 40ºC

Detection : UV 280 nm

Injection Vol. : 1 µL

Sample : Fat-soluble Vitamins

Solutes :1. Vitamin K3 (0.14 mg/mL, tR 3.703)
2. Vitamin A (0.14 mg/mL, tR 5.033)
3. Vitamin A acetate (0.14 mg/mL, tR 5.740)
4. Vitamin D3 (0.14 mg/mL, tR 6.943)
5. Vitamin E (0.14 mg/mL, tR 9.027)
6. Vitamin K1 (0.14 mg/mL, tR 9.750)