
Vitamins



Column : Inertsil ODS-3V
5µm
150 × 4.6 mm I.D.

Column Cat. No. :

Eluent : CH_3OH

Flow rate : 1.0 mL/min.

Col. Temp. : 40ºC

Detection : UV 280 nm

Injection Vol. : 1 µL

Sample : Fat-soluble Vitamins

Solutes : 1. Vitamin K3 (0.14 mg/mL, tR 2.130)
2. Vitamin A (0.14 mg/mL, tR 3.103)
3. Vitamin A acetate (0.14 mg/mL, tR 4.227)
4. Vitamin D3 (0.14 mg/mL, tR 6.650)
5. Vitamin E (0.14 mg/mL, tR 8.063)
6. Vitamin E acetate (0.14 mg/mL, tR 10.820)
7. Vitamin K1 (0.14 mg/mL, tR 15.270)