
Vitamins



Column : Inertsil CN-3
5µm
250 × 4.6 mm I.D.

Column Cat. No. :

Eluent : n-Hexane/Ethanol = 99.5/0.5

Flow rate : 1.0 mL/min.

Col. Temp. : 40ºC

Detection : 0 - 3.3 min: VIS 450 nm,
3. 3 - 17 min: UV 280 nm

Injection Vol. : 3 µL

Sample : Fat-soluble Vitamins

Solutes :1. β - Carotene (0.13 mg/mL, tR 2.890)
2. Vitamin K1 (0.07 mg/mL, tR 3.557)
3. Vitamin A acetate (0.07 mg/mL, tR 3.833)
4. Vitamin K3 (0.07 mg/mL, tR 5.843)
5. Vitamin E (0.07 mg/mL, tR 8.047)
6. Vitamin D3 (0.07 mg/mL, tR 10.193)
7. Vitamin A (0.07 mg/mL, tR 14.927)