
Insulin standards E10



Column	: Inertsil ODS-2 5µm 250 × 4.6 mmI.D.
Column Cat. No.	: 5020-01102
Eluent	: A:0.1% CF_3COOH in CH_3CN B:0.1% CF_3COOH A/B 20/80?(5min)?30/70 →(40min)→50/50
Flow rate	: 1.0 mL/min.
Col. Temp.	: RT
Detection	: UV 280 nm
Injection Vol.	:
Sample	:
Solutes	:porcine insulin Bovine insulin