
Purification of Zingiberis rhizoma



Column	: Inertsil PREP-ODS 10µm	Solutes	: 1. Zingiberis rhizona
Column Cat. No.	: 5020-35151		
Eluent	: A:CH_3OH B:H_2O A/B 60/40→(35min)→100/0		
Flow rate	:		
Col. Temp.	:		
Detection	: UV 254 nm		
Injection Vol.	:		
Sample	: Zingiberis rhizona 250 × 6.0 mmI.D. - Column size - 250 × 50 mmI.D. 1mL/min - Flow rate - 70mL/min 0.5mL - Sample size - 35mL		