
Analysis of Betamethasone valerate



Rs=14.2

Column : Inertsil ODS-3V
5µm
250 × 4.6mmI.D.

Column Cat. No. :

Eluent : CH_3OH/H_2O=70/30

Flow rate : 1.0 mL/min.

Col. Temp. : 25°C

Detection : UV 254 nm

Injection Vol. : 10µL

Sample : 1. β - methasone zalerate(0.05mg/mL)
2.Isoamyl benzoate(0.5µg/mL)

Solutes : 1. Betamethasone valerate (0.05 mg/mL)
2. Isoamyl benzoate(0.5µg/mL)