
Analysis of Idoxuridine



Column	: Inertsil ODS-3V 5µm 150 × 4.6mmI.D.
Column Cat. No.	:
Eluent	: CH_3OH/H_2O=4/96
Flow rate	: 1. 3 mL/min.
Col. Temp.	: 30ºC
Detection	: UV 254 nm
Injection Vol.	: 10µL
Sample	: 1. 2'-Deoxyuridine (4ng/mL) 2. Iodouracil(12ng/mL) 3. Idoxuridine (20ng/mL)

Solutes	: 1. 2'-Deoxyuridine (4ng/mL) 2. Iodouracil(12ng/mL) 3. Idoxuridine (20ng/mL)