


Tocopherol

1. Tocopherol (1.05µg)



2. Tocopherol acetate (1.20µg)





Column : Inertsil ODS-3
5µm
150 × 4.6mmI.D.

Column Cat. No. :

Eluent : MeOH / H_2O (98 : 2)

Flow rate : 1.0 mL/min.

Col. Temp. : 30°C

Detection : UV 292 nm

Injection Vol. :

Sample : 1. Tocopherol(1.05µg)
2. Tocopherol acetate(1.20µg)

Solutes : 1. Tocopherol(1.05µg)
2. Tocopherol acetate(1.20µg)