


Procaine







Column : Inertsil ODS-3
5µm
150 × 4.6mmI.D.

Column Cat. No. :

Eluent : 50mM KH2PO4 (pH3.0) / MeOH / Sodium 1-pentanesulfonate
(800 : 200 : 1)

Flow rate : 1.0 mL/min.

Col. Temp. : 40ºC

Detection : UV 254 nm

Injection Vol. :

Sample : 1. Procaine Hydrochloride (1.25g)
2. Caffeine (1.25 g)

Solutes : 1. Procaine Hydrochloride (1.25g)
2. Caffeine (1.25 g)