
Estrogens



Column : Inertsil ODS
5µm
250 × 4.6 mmI.D.
Column Cat. No. : 5020-02202
Eluent : CH_3CN/1% $CH_2ClCOOH$(pH3 by NH_4OH)=60/40
Flow rate : 1.0 mL/min.
Col. Temp. : 40°C
Detection : UV 254 nm
Injection Vol. :
Sample :

Solutes : 1. Estriol 4mg/mL
2. Estradiol 4mg/mL
3. Estrone 2mg/mL