
Catechol amine



Column : Inertsil C4
5µm
150 × 4.6mmI.D.

Column Cat. No. : 5020-01401

Eluent : CH_3CN/5mM IPCC-12A+20mM NaH_2PO_4(pH3 by H_3PO_4)=20/80
IPCC-12A:Sodium Dodecylsulfate

Flow rate : 1.0 mL/min.

Col. Temp. : 40ºC

Detection : UV 280 nm

Injection Vol. :

Sample :

Solutes : 1. Vanillylmandelic acid(VMA)
2. 3,4-Dihydroxy phenyl acetic acid
3. Homovanillic acid(HVA)
4. L-DOPA
5. L-Noradrenaline
6. L-Adrenaline
7. Dopamine Hydrochloride