
Analysis of Adrenocortical hormones



Column	: Inertsil ODS-3V
	5µm
	150 × 4.6mmI.D.
Column Cat. No.	:
Eluent	: Methanol/Water=60/40
Flow rate	: 1.0mL/min
Col. Temp.	: 40ºC
Detection	: UV 254 nm
Injection Vol.	: 3µL
Sample	: Prednisone (0.1mg/mL)
	Cortisone (0.1mg/mL)
	Prednisolone (0.1mg/mL)
	Corticosterone (0.1mg/mL)

Solutes : 1. Prednisone
2. Cortisone
3. Prednisolone
4. Corticosterone