
# Analysis of Estrogens




| | |
|---|---|
| Column | : Inertsil CN-3<br>5µm<br>250×4.6mmI.D. |
| Column Cat. No. | : |
| Eluent | : Acetonitrile / Water = 35 / 65 |
| Flow rate | : 1.0mL/min |
| Col. Temp. | : 25ºC |
| Detection | : UV 220 nm |
| Injection Vol. | : 1µL |
| Sample | : Estriol (1mg/mL)<br>β - Estradiol (1mg/mL)<br>Estrone (1mg/mL)<br>Ethynylestradiol (1mg/mL)<br>Hexestrol (0.5mg/mL)<br>Diethylstilbestrol (0.5mg/mL) |
| Solutes | : 1. Estriol<br>2. Estradiol<br>3. Estrone<br>4. Ethynylestradiol<br>5. Hexestrol<br>6. Diethylstilbestrol |