
Pesticides



Column : Inertsil C4
5µm
150 × 4.6mmI.D.

Column Cat. No. : 5020-01401

Eluent : THF/10mM KH_2PO_4+0.1% EDTA-2Na
(pH3.5 by H_3PO_4)=70/30

Flow rate : 1.0 mL/min.

Col. Temp. :

Detection : UV 254 nm

Injection Vol. :

Sample :

Solutes : 1. Asulam 130µg/mL
2. Oxine-copper 40µg/mL
3. Acetophenone 0.14µg/mL
4. Thiuram 190µg/mL