
Analysis of Carbamates



No.	R.Time	Peak Name	Area
1	2.080	Oxamyl	7.08845e+06
2	2.323	Methomyl	5.65873e+06
3	3.944	Aldicarb	2.97407e+06
4	4.945	Bendiocarb	7.58038e+06
5	6.456	Ethiofencarb	1.3242e+07
6	7.312	Pirimicarb	1.29822e+06
7	11.568	Fenobucarb	6.53701e+06
8	13.359	Methiocarb	7.07763e+06

Column : Inertsil ODS-3
5µm
150 × 4.6mmI.D.

Column Cat. No. :

Eluent : 53% MeOH

Flow rate : 1.0 mL/min

Col. Temp. : 40ºC

Detection : FL Ex 345nm Em 450nm
POST REACTION NaOH 100? to OPA 50?

Injection Vol. : 25µL

Sample : Pesticides

Solutes : 1 . Oxamyl 2.83mg
2 . Methomyl 2.5mg
3 . Aldicarb 1.65mg
4 . Bendiocarb 2.5mg
5 . Ethiofencarb 5.10mg
6 . Pirimicarb 7.8mg
7 . Fenobucarb 2.0mg
8 . Methiocarb 2.4mg
in 50mL 52% MeOH