
Analysis of Phthalic acid dialkylester



No.	R.Time	Area
1	2.606	213991
2	3.162	159440
3	4.173	288422
4	6.017	75183
5	9.798	97919
6	11.911	84698
7	15.092	70777
8	17.003	58989

Column : Inertsil C8-3
(5µm , 150 × 4.6 mmI.D.)
Column Cat. No. : 5020-01900
Flow Rate : 1.0mL/min
Col. Temp. : 40ºC
Detection : UV 280nm
Eluent : A:H O
B :CH_3CN
A/B=25/75 -(5min)- 25/75 -(10min)- 0/100,v/v
Sample : Phthalic acid dialkyl esters

Solutes : 1. Diethyl phthalate(1 µL/mL)
2. Di-iso-propyl Phthalate(1µL/mL)
3. Di-n-butyl Phthalate(1 µL/mL)
4. Di-iso-butyl Phthalate(1 µL/mL)
5. Dicyclohexyl Phthalate(0.5 mg/mL)
6. Di-n-heptyl Phthalate(1 µL/mL)
7. Di-(2-ethylhexyl) Phthalate(1 µL/mL)
8. Di-n-decyl Phthalate(1 µL/mL)
9. Di-n-dodecyl Phthalate(1 µL/mL)