
Analysis of Amitriptyline



Column	: Inertsil C8-3
	150 × 4.6 mmI.D.
Column Cat. No.	: 5020-01900
Eluent	: Methanol/20mM Phosphate buffer(pH7.0)=65/35
Flow rate	: 1.0 mL/min.
Col. Temp.	: 40ºC
Detection	: UV 254 nm
Injection Vol.	: 5µL
Sample	:

Solutes : 1. Uracil(0.014mg/mL)
2. n-Butyl Paraben(0.020mg/mL)
3. Propranolol(0.48mg/mL)
4. Dipropyl Phthalate(0.20mg/mL) + Naphthalene(0.10mg/mL)
5. Acenaphthene(0.19mg/mL)
6. Amitriptyline(0.09mg/mL)