
Analysis of Polystyrene Standard



Column	: Inertsil SIL-100A
	5 µm
	250 x 4.6 mm I.D.
Column Cat. No.	: 5020-01712
Eluent	: THF
Flow rate	: 1.0 mL/min.
Col. Temp.	: Ambient
Detection	: UV 254 nm
Injection Vol.	: 5 µL
Sample	:

Solutes : 1. Polystyrene (0.2 mg/mL, Mp 1,020,000)
2. Polystyrene (0.2 mg/mL, Mp 115,000)
3. Polystyrene (0.2 mg/mL, Mp 9,200)
4. Polystyrene (0.2 mg/mL, Mp 950)