
Analysis of Planar and Non-planar Compounds



Column	: Inertsil ODS-3	Solutes	: 1. Benzo[a]pyrene (BaP, Planar)
	5 µm		2. Phenanthro[3,4-c]phenanthrene (PhPh, Non-planar)
	250 x 4.6 mm I.D.		3. Tetrabenzonaphthalene (TBN, Non-planar)
Column Cat. No.	: 5020-01732		
Eluent	: A) CH_3CN		
	B) H_2O		
	A/B = 85/15, v/v		
Flow rate	: 2.0 mL/min.		
Col. Temp.	: 30°C		
Detection	: UV 254 nm		
Injection Vol.	:		
Sample	:		