
Analysis of sweeteners



System : GL-7400 system

Column : Inertsil NH2 (5µm, 150×4.6 mm I.D.)

Column Cat. No. :

Eluent : A) CH_3CN
B) H_2O
A/B = 80/20, v/v

Flow rate : 1.0 mL/min

Col. Temp. : 40 °C

Detection : UV 210 nm (GL-7452A PDA detector)

Injection Vol. : 20 µL

Sample : Sweeteners

Solutes : 1. Stevioside(10 mg/L)
2. Rebaudioside A(10 mg/L)