
Hydrocarbons C2-C4

Specification
One GC Channel



Conditions

System	**:** Micro GC
Column	**:** PoraPLOT Q 10 m
Column Temp.	**:** 100 ℃
Carrier Gas	**:** He 190 kPa
Inject Time	**:** 40 msec

Analyte:

1. Ethylene (C_2H_4)	1000 ppm (v/v)
2. Ethane (C_2H_6)	1000 ppm (v/v)
3. Propylene (C_3H_6)	1000 ppm (v/v)
4. Propane (C_3H_8)	1000 ppm (v/v)
5. Isobutane (C_4H_{10})	1000 ppm (v/v)
6. *n*-Butane (C_4H_{10})	1000 ppm (v/v)

Balance Gas N_2