
Hydrocarbons C1-C2

Specification
One GC Channel



Conditions

System	**:** Micro GC
Column	**:** PoraPLOT U 10 m
Column Temp.	**:** 40 ℃
Carrier Gas	**:** He 200 kPa
Inject Time	**:** 40 msec

Analyte:

1. Methane (CH_4)	100 ppm (v/v)
2. Carbon dioxide (CO_2)	
3. Ethylene (C_2H_4)	100 ppm (v/v)
4. Ethane (C_2H_6)	100 ppm (v/v)
5. Acetylene (C_2H_2)	100 ppm (v/v)

Balance Gas N_2